

HemeRx™

HEMOSTATIC COTTON PELLETS WITH EPINEPHRINE

HemeRx pellets are cotton pellets saturated with a solution of racemic epinephrine hydrochloride and dried using a technique that retains the efficacy of the medication. Normal tissue moisture activates the medication, promptly retracting gingival tissue and offering local hemostasis without any tissue damage.

Indications for Use: For use in dental practice whenever gingival retraction is advisable, e.g. preceding impression taking with any elastic impression material, prior to placement of caps, crowns, bridges or other restorative materials, or during tooth preparation. HemeRx pellets will help control capillary oozing and tissue engorgement through temporary vasoconstriction. Normal gingival moisture will activate HemeRx by facilitating the absorption of the medication. Otherwise, the pellets may be quickly dipped in water just before being applied.

Contraindications: HemeRx pellets are contraindicated for persons who are epinephrine sensitive or for persons suffering from cardiovascular disease, diabetes, or hyperthyroidism. Caution should also be taken not to exceed the recommended amount for medically compromised, elderly, pregnant, and pediatric patients.

Warnings: Contains Sodium Bisulfite, a sulfite that may cause allergic-type reactions, including anaphylactic symptoms and life-threatening or less severe asthmatic episodes in certain susceptible persons. The overall prevalence of sulfite sensitivity in the general population is unknown, but probably low. Sulfite sensitivity is seen more in asthmatic than in non-asthmatic persons.

Precautions:

Read and follow all instructions.

For professional in-office use only.

For topical intraoral use only. Do not ingest.

Keep out of reach of children.

Do not use on patients with epinephrine allergy or sensitivity.



Manufactured for

Sprig Oral Health Technologies

6140 Horseshoe Bar Road, Suite L

Loomis, CA 95650 USA

Tel: +1 916.677.1447 | www.sprigusa.com



SN70301/0417